Diamond Shield IE

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Mikroströme - EMS

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Operation.....41

Technical Data.....48

Manufacturer’s warranty declaration.....49

For disposal.....49

Accessories.....50

Declaration of Conformity.....53

Introduction

The Diamond Shield IE Zapper

The Wellness Device with Microcurrents


- wide frequency range of 0.5 Hz to 2.5 MHz.
- intensity range from small voltages of 0.1 V up to high voltages of 16 V
- Impulse discharge
- Wobble functionality
- Modulation of the output frequencies
- Constant current control
- Running frequencies from special chip cards (optional accessory)

The operation can be done via a 9V battery or a rechargeable battery.

The device can be used to enhance the well-being.

An earthing cable and an enclosed power socket tester are included in the scope of delivery.

The good effect of the earthing (so called negative ions) occurs only if the power socket is properly earthed. Be sure to check this with the power socket tester.

Warning

The use of the earthing plug (or any other electrical appliances) to an incorrect power socket can be perilous.

Do not use the earthing cable during a thunderstorm.

Important Notices

Safety Instruction

- Use only the supplied grounding cable. Under no circumstances use a different cable for grounding.
- Under any circumstances never insert a plug into an electrical outlet.
- Always insert first the grounding cable to the Diamond Shield and then connect it to the power outlet.
- Never use this device to anyone who has a pacemaker.
- Never use this device in the eye area.
- Keep the Diamond Shield always clean and dry – wiping it with a damp cloth is permitted.

***Using of the device according to the operating instructions is essential.
Mannayan GmbH & Co KG assumes no liability for misuse.***

Note

In case of feeling a tingle or burning on your skin during the usage please lower the intensity (using the arrow down key). The efficiency of the application is independent from the above. When hypersensitive to metal, we gladly exchange the metal bracelets with cloth bracelets.

Unpacking

Before using this device, please check the completeness of the supplied accessories:

- Connection cable for bracelets or other electrodes (incl. two pairs of connector adapters)
- 2 metal bracelet electrodes
- 2 adhesive electrodes
- Grounding cable and socket tester
- 2 Chipcards
- 9V Battery

Storage and transportation

The storage of the equipment shall be always stored in dry, closed rooms. If the device is transported or stored in extreme temperatures, a period of at least 2 hours should be observed for acclimatization before handling.

The device is intended for use in clean, dry rooms. They cannot be used in extreme dust, damp, dusty air, danger of explosion or chemical vapors. Protect the Diamond Shield IE from direct sunlight.

The permissible operating temperature ranges from +15°C to +30°C. During storage or transportation, the ambient temperature can be between -20°C and +70°C. If condensation has been formed during transportation or storage, the device must acclimate for about 2 hours and dried by appropriate circulation. Thereafter, the operation of the device is permitted. Nominal data with tolerances are valid after a warm-up time of at least 30 minutes at an ambient temperature of 23°C. A value without tolerances is a guiding value for the device.

Maintenance


The outside of the device should be cleaned regularly with a soft, non-linting cloth. The Diamond Shield IE should never be cleaned with alcohol or other solvents.

Intended Usage

1. The device can be used to improve the subjective well-being of people who are open to such an approach.
2. It is applicable as EMS-device.
3. With a connected speaker it can serve for acoustic frequency application (music therapy).
4. The preinstalled Diamond Shield program can be used to balance the chinese meridian system.

Feature Identification

Front Side

- | | | |
|----------------|------------------|---------------|
| 1. LED | 4. START Button | 7. OK Button |
| 2. LCD Display | 5. On-Off Switch | 8. Arrow Keys |
| 3. Mute Switch | 6. EXIT Button | |

Connection Side

1. Grounding connection
2. Chipcard slot
3. Connection for electrodes (ie bracelets)

Back Side

1. Battery compartment

Bringing into service

Insert the battery


1. Open the battery compartment cover on the back of the device.
2. Insert the supplied battery. Please keep in mind that the battery has to be inserted correctly. Meaning the plus has to go to plus and minus to minus. (If however the battery is inserted incorrectly, the device will not be damaged by that.)

We recommend using alkaline batteries or rechargeable batteries.

Connecting the cables

Connect the electrode cable for the bracelet electrodes in the right socket (above LED). Then connect the bracelets.

If you want to work with grounding: Connect the grounding cable to the left socket. Then connect the ground wire to the outlet.

Operation

Switch On by pressing briefly the **POWER** button

Switch Off by pressing the **POWER** button and holding it for about 1 second

Automatic Switch Off

After 12 minutes without user or program activity the device switches off automatically to save battery power.

Charge condition of the battery

The state of charge of the battery is displayed on the LED.

If the device is switched on, the LED flashes. When the battery is full, the flashing frequency is slow. When the battery is almost empty the flashing frequency is faster.

If the battery is completely empty, a message appears, and the device can no longer be used.

Activate / deactivate the loudspeaker

Press the **SOUND** button to switch the loudspeaker on or off. The current state is then displayed briefly.

Keypad lock


The device has a keypad lock to prevent accidental operation of buttons. To activate or deactivate the keypad lock briefly press the **POWER**-button and then press **OK**. The device will display the status of the keypad lock then.

Menu Navigation

Navigation through the menu structure of the Diamond Shield IE is simple and identical for all functional areas:

	The currently selected item is marked with a small black rectangle. A scroll bar on the right edge display appears, when the current list has more entries than can be displayed. The scroll bar then indicates, where you are in the list.
--	---

Use the **up-arrow-key** **+** and the **down-arrow-key** **-** to navigate to the desired item.

For switching quickly between the first and the last entry you can jump from the last to the first item with the **down-arrow-key**, and accordingly from the first to the last item with the **up-arrow-key**.

Confirm your selection with the **OK** button. Then a submenu appears that corresponds to the selected item.


Press the **EXIT** key to return one level up (if you are not already in the main menu).

Important Terms

The structure of the application programs from the standard menu or from chip-cards is composed of two levels of elements:

- **Frequency program:** A frequency program is a defined frequency application with all its properties such as the frequency, wobble, intensity, modulation and duration.
- **Sequence:** A sequence is a succession of several single frequency programs. For example, there is a sequence with the name „Diamond Shield“, which includes the frequency programs „DS1“, „DS2“, ... up to „DS26“.


Menu „Main menu“


Shortly after switching on the device the **Main menu** will appear. You can always return to the **Main menu**, if you want to leave the submenus, by pressing the **EXIT** button.

The Main menu shows the following items:

- **Chipcard:** This item is displayed only if a chip card is inserted. It can be used to display and run the contained programs.
- **Standard:** Under this menu item you can find preinstalled sequences and frequency programs.
- **Settings:** In this menu various settings (language, display contrast, constant amps, speaker) can be undertaken. (see page 45)


„Standard“ menu

Under the menu item **Standard** you will find the pre-installed frequency programs, with which the Diamond Shield IE Zapper is supplied by default. Most of these programs consist of a sequence of several sub-programs.

The first level of the menu shows the available sequences. When a sequence is selected and then the **OK** button is pressed, the frequency programs of the selected sequence are shown.


Executing a sequence

Select the desired sequence with the **+** and **-** button. To run the entire sequence press the **START** button.


First the device will test the electrical contact (see the chapter „Contact control“). Then all the individual programs of the selected sequence are successively executed.

Running a single program of a sequence

Select the desired sequence with the **+** and **-** button. Press the **OK**-button to enter the submenu, where the sub-programs of this sequence are displayed.

Then select the desired subprogram and run it by pressing the **START**-button. Previously the device checks the electrical contact (see **Contact control**).

Contact control

Before a sequence or a single program is executed, the device checks whether the required electrical contact is made. If the contact is not good enough, you will hear a beep every 3 seconds. By pressing the **OK**-button, the frequency program can be executed anyway. However no effect will occur.

If the required contact is made, the display shows for a short time the resistance in k-Ohms and then automatically begins executing the frequency program.

Display

Diam.Shield	DS 1
6:02	5:56
0.6V 0.30mA	0.2V

- In the headline the sequence name and the name of the current single program are shown in turns.
- Below is the remaining execution time displayed (minutes: seconds).
- Below the intensity is displayed.
 - For frequency programs with constant current control, the actual voltage is displayed on the left side. On the right side the selected intensity is displayed in milliamps (1/1000 amps).
 - For frequency programs without current control, only the voltage in volts is displayed.

Operating options

Pressing the **START**-button will interrupt the program.

Repressing the **START**-button will play the program from the point of interruption.

The intensity can be adjusted using the arrow keys **up** and **down**.


Pressing the **EXIT**-button will quit the current application.

Running a chip card

Chip cards are available as accessories containing preprogrammed frequency sequences. Both these diamond chip cards as well as individual chip cards can be played as follows:

Insert the chipcard into the chipcard slot of the device.
Select the menu item **Chip card** and press the **OK**-button.


A submenu appears, which displays the contents of the chip card. The entire sequence(s) and/or sub-programs of **Chip card** can be applied similarly as it is described in the chapter „**Standard“ menu**.

For running the entire sequence of the selected program, press the **START**-button


For running single sub-programs of the sequence, enter the submenu of the sequence by pressing the **OK**-button.

The subprograms of this sequence are displayed. Then select the desired subprogram and run it by pressing the **START**-button.

Menu „Config“

This menu is used to configure the device. The settings made here are retained even after switching off the device or changing the battery.

- **Language:** Select from English or German.
The input is done with the arrow keys **left** and **right**.
- **Contrast:** Adjusts the contrast of the LCD display.
The input is done with the arrow keys **left** and **right**.
- **Const amps:** enables or disables constant current control.
The input is done with the arrow keys **left** and **right**.
If **no** is selected, the intensity of all the frequency programs (including those, that are stored on the chip cards) is usually set on the voltage and does not have a constant current control.
- **Speaker:** the tone pitch of the internal speaker.
Use the arrow keys **left** and **right** to choose 4 or 5 kHz.
- **HW Rev.:** Displays the hardware version of your device.

Technical Data

General

- Operating temperature: +15° C to +30° C
- Storage temperature: -20° C to + 70° C
- Power supply: 9V alkaline (manganese) battery or comparable accumulator.

Connections

- 3.5 mm socket (stereo) for grounding plug
- 3.5 mm socket (stereo) for signal outputs
- Chip card slot

Characteristic of the output signals

- Frequency range 0.5 Hz - 2,5 MHz
- Waveform: square wave signal
- Minimum accuracy of the frequencies: 100 ppm
- Output Voltage: 0.1 - 16 Volts (measured from the 0-Volt-line)
- Offset: No offset, positive offset, positive offset + 0.3 Volts
- Modulation up to 254 times the output frequency (but 2,5 MHz must not be exceeded)
- Wobble (maximum up to twice the fundamental frequency)
- Voltage control to keep the output current constant (only at frequencies above 20Hz)
- Internal resistance: 1 kOhm

Manufacturer's warranty declaration

Mannayan GmbH & Co KG grants a two-year manufacturer's warranty for the purchased product from the date of purchase. Warranty claims are valid throughout all of Europe. There is no restriction of warranty claims as a result of the manufacturer's warranty conditions. This guarantees a right to the coverage of repair and shipping costs, which occur as a result of a defect during the warranty period.

However, the warranty claim only remains valid if the defect is not due to improper handling or normal wear and tear. In the event of a warranty claim, you must contact the named company. You will then receive documents which will enable you to return the goods free of charge. Should the warranty claim arise outside Europe, the device must be returned to the manufacturer at the customer's expense. If a repair is no longer possible, the manufacturer will exchange it for an equivalent product.

For disposal

The Diamond Shield IE device must – as well as other electronic devices – not be disposed with normal household waste. You can bring it to the local separate collection system for electrical appliances or return it to the address below (unfree postage will not be accepted)

Please also note, that the disposal of batteries is prohibited on the normal household waste for environmental reasons. You can bring them to the public collection places or wherever batteries are sold.

You are welcome to return the battery to the following address (unfree postage will not be accepted):

Mannayan GmbH & Co. KG
Unterer Anger 15/Rgb.
80331 Munich
Germany


Accessories

Included in the scope of delivery:


Double cable to connect the electrodes to the Zapper

The enclosed adapters (10mm) are compatible to all electrodes excepting different adhesive and textile wristband electrodes. The cable with the red plug is the current carrying line (positive pole). For plate zapping, the plate is interconnected in this line.

Metal Bracelet Electrodes

The width of the metal bracelet electrodes can be adjusted. For this purpose the blue element can be opened at the indicated edge.

Adhesive electrodes

The adhesive electrodes can be placed on the body for specific local application.

Earthing cable

The insert pins of the earthing cable plug are non-conductive. Only the two lateral earthing contacts are conductive. The cable must first be connected to the zapper, afterwards to a grounded outlet. The grounding capacity of the socket can be checked by means of the supplied socket tester.

During the course of the program the earthing occurs via the impulse discharge function. It is recommended to stay earthed for about 50 minutes after completion of the program. For this purpose the zapper can be switched off.

Optional:


Fractal Enhancer with more than 1200 fractals, arranged as Tesla coils. It can be interconnected within the current carrying line by adding a connecting cable.

Wristband Electrodes with sponge

The adjustable wristbands need to be moistened for use.

Adhesive Electrodes, small (45 x 35 mm) for small area applications

9-Volt-Akku

a 9-Volt-Lithium-Akku with a suitable charging cable



Hand Electrodes, stainless steel

The hand electrodes are versatile and easy to clean.

Textile electrodes with silver fibre are available as wrist cuffs, elbow cuffs, knee cuffs, socks and gloves. They are soft, breathable and feel comfortable on the skin.

Optional:

Surface or fabric electrodes are made of current-conducting textile and can be used for large-area application on humans and animals. They are equipped with a velcro fastener or a fastening strap.

Chipcards

A constantly growing range of available chip cards complements the zapper and expands its functionality to a great extent.

Leather etui

Genuine leather etui with 32 compartments for sorting and storing your chipcards. The etui fits into the zapper case.

Plate, stainless steel

For plate zapping, the plate is interconnected within the current carrying line (the cable with the red plug).

Zappicator

The zappicator is connected to the zapper with the current carrying cable. Then you can place food on the plate. In the Standard menu you can find the zappicator program. The earthing cable is not needed.

Speaker

The speaker can be used for acoustic frequency application.

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© June 2024 Mannayan GmbH & Co. KG
All rights reserved
Any kind of reproduction only with
approval of Mannayan GmbH & Co. KG

Mannayan GmbH & Co. KG
Dorfstraße 4a
84387 Julbach
Germany

Tel: +49 (0) 89-18 90 84 6-0
Fax: +49 (0) 89-18 90 84 6-15
E-Mail: info@mannayan.com
Web: <http://www.mannayan.com>

Self-made chip cards for the most popular Zapper models in just a few minutes!

The screenshot shows the 'Chipcard Creator' software interface. It features a 'Library' tab and an 'Edit programs' button. A search bar is present with the text 'Select the libraries to search:'. Below this is a table with columns for 'Library', 'Indications', and 'Also for'. The table lists various chipcards and their associated conditions and frequencies.

Library	Indications	Also for
Chipcard	Veneral disease according to Dr. Clark	AIDS HIV
Chipcard	ALS Amyotrophic Lateral Sclerosis according to Dr. Clark	ALS Amyotrophic Lateral Sclerosis
Chipcard	5 Elements Fire Support according to Baklayan	Abdomen - pain

The 'Also for' column lists various conditions such as Veneral disease, Syphilis, Treponema pallidum, Wart, Human Papilloma, Wart HPV, Cervical smear, Herpes simplex 1, Herpes simplex 2, Herpes simplex 3, Hepatitis B, Chlamydia, Neisseria gonorrhea, Trichomonas, Haemophilus influenzae, Fasciola hepatica, Clonorchis sinensis, Fasciolopsis buski, Eurytema pancreaticum, Shigella 1, Shigella 2, Shigella 3, Nocardia 1, Nocardia 2, Echinostoma revolutum, Echinoparyphium recurvatum, Gastrothylax elongatus, Streptococcus 1, Streptococcus 2, Streptococcus 3, Staphylococcus, 5 Elements Fire Support, Concentration power reduced, Diarrhea, Insomnia, Palpitation, Perspiration - spontaneous, Reduced ability to remember, Stool with undigested residues, and Precordial anxiety.

Overlaid on the screenshot is a white USB Zapper device with a blue circular logo. Three white chipcards are also shown, each labeled 'Healing Frequency' and 'Diamond Shield'.

www.healingfrequency.com

THE NEW ONLINE LIBRARY FOR HARMONIC AND HEALING FREQUENCIES

Create your own, individual smart cards easily, fast and conveniently on your home computer!